



Beagle Lateral bedmover



*Lift your products to the
next level of controlling*



Key Features

- 360-degree movement: all directions motor
- drive system
- Highly efficient and safe electronic controlling system
- Foldable steering device (optional)
- Extreme compact connecting possible for best fitting in elevators
- Three-phase highly efficient and very powerful flange motor inside
- Steering possible at small patient rooms
- Large battery capacity
- Sophisticated customized software settings available

Technology

Advanced processor technology for smooth steering and personal settings result in safety during driving for both patient and driver.

Dimensions:

790mm (L) x 550mm (W) x 390mm-1120mm (H)

Speed: up to 10km/h (adjustable)

Load capacity: 1000 kg

Batteries: 80 Ah 24 volt

Charger: 8 Amp 24 Volt IP68

Weight: 160 kg

Beagle Mobility BV

Oenerweg 20
8161 PM Epe, The Netherlands
Phone : +31 578 690 457
E-mail : info@beagle.nl
Website : www.beagle.nl

VAT : NL 823143193B01
Chamber of commerce : 51196875

Swift : RABONL2U
IBAN : NL33rabo0129267287